
Germline SNV and Indel Calling

10/26/20


Sam Peters and Chris Miller

How do we identify Germline SNVs and Indels?

- 1. Align Reads to Reference
- 2. Call Genotypes
- 3. Annotate variants
- 4. Filter Variants
 - 1. Annotations
 - Mapping Quality (MQ)
 - 2. Read Depth (DP)
 - 3. Genotype Quality (GQ)
 - 2. Allele Frequency
 - 3. Region of Interest
- 5. Final QC
 - 1. Manually Check Sequence Context (IGV)
 - Homopolymers, Repeat Regions are hard to sequence




Start with Analysis-Ready Reads




Read Alignment Visualization

Exome Sequencing Data from NA12878



Genotype NGS samples with GATK HaplotypeCaller



Call Genotypes Using GATK HaplotypeCaller

How HaplotypeCaller works

1. Define active regions

The program determines which regions of the genome it needs to operate on (active regions), based on the presence of evidence for variation.

2. Determine haplotypes by assembly of the active region


For each active region, the program builds a De Bruijn-like graph to reassemble the active region and identifies what are the possible haplotypes present in the data. The program then realigns each haplotype against the reference haplotype using the Smith-Waterman algorithm in order to identify potentially variant sites.

3. Determine likelihoods of the haplotypes given the read data

For each active region, the program performs a pairwise alignment of each read against each haplotype using the PairHMM algorithm. This produces a matrix of likelihoods of haplotypes given the read data. These likelihoods are then marginalized to obtain the likelihoods of alleles for each potentially variant site given the read data.

4. Assign sample genotypes


For each potentially variant site, the program applies Bayes' rule, using the likelihoods of alleles given the read data to calculate the likelihoods of each genotype per sample given the read data observed for that sample. The most likely genotype is then assigned to the sample.



https://gatk.broadinstitute.org/hc/en-us/articles/360035531412

https://gatk.broadinstitute.org/hc/en-us/articles/360037225632-HaplotypeCaller

Two Methods for Variant Calling with HaplotypeCaller



Unfiltered VCF Output

header

```
#contig=<ID=HLA-DRB1*01:02:01,length=11229>
#contig=<ID=HLA-DRB1*03:01:01:01,length=13908>
##contig=<ID=HLA-DRB1*03:01:01:02,length=13426>
#contig=<ID=HLA-DRB1*04:03:01,length=15246>
#contig=<ID=HLA-DRB1*07:01:01:01,length=16110>
##contig=<ID=HLA-DRB1*07:01:01:02,length=16120>
#contig=<ID=HLA-DRB1*08:03:02,length=13562>
#contig=<ID=HLA-DRB1*09:21,length=16039>
#contig=<ID=HLA-DRB1*10:01:01,length=13501>
#contig=<ID=HLA-DRB1*11:01:01,length=13921>
#contig=<ID=HLA-DRB1*11:01:02.length=13931>
#contig=<ID=HLA-DRB1*11:04:01,length=13919>
#contig=<ID=HLA-DRB1*12:01:01,length=13404>
#contig=<ID=HLA-DRB1*12:17,length=11260>
#contig=<ID=HLA-DRB1*13:01:01,length=13935>
#contig=<ID=HLA-DRB1*13:02:01,length=13941>
##contig=<ID=HLA-DRB1*14:05:01,length=13933>
#contig=<ID=HLA-DRB1*14:54:01,length=13936>
#contig=<ID=HLA-DRB1*15:01:01:01,length=11080>
#contig=<ID=HLA-DRB1*15:01:01:02,length=11571>
#contig=<ID=HLA-DRB1*15:01:01:03,length=11056>
#contig=<ID=HLA-DRB1*15:01:01:04,length=11056>
 contig=<ID=HLA-DRB1*15:02:01,length=10313>
#contig=<ID=HLA-DRB1*15:03:01:01,length=11567>
#contig=<ID=HLA-DRB1*15:03:01:02,length=11569>
##contig=<ID=HLA-DRB1*16:02:01,length=11005>
#source=HaplotypeCaller
#bcftools_viewVersion=1.10.2-91-g365d117+htslib-1.10.2-109-gdcd4b73
#bcftools_viewCommand=view -r chr13 NA12878-HG001-merged.vcf.gz; Date=Tue Oct 20 08:55:27 2020
```

calls


```
##bcftools_viewCommand=view -r chr13 NA12878-HG001-merged.vcf.gz; Date=Tue Oct 20 08:55:27 2020
#CHROM
       POS
                                                 FILTER INFO
                                                                  FORMAT NA12878-HG001
                                 ALT
chr13
        16002338
                                                 61.65
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=-0.674; DP=4; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=25.51; MQRankSum=0.319; QD=15.41; ReadPosRankSum=-0.319; SQR
        16003734
                                                 64.64
                                                                  AC=1; AF=0.5; AN=2; BaseORankSum=-0.842; DP=6; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MO=48.15; MORankSum=0.842; DD=10.77; ReadPosRankSum=-0.842; SOR
chr13
                                         G
                                 G
        16003747
                                                 61.64
                                                                  AC=1; AF=0.5; AN=2; BaseORankSum=-0.366; DP=7; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MO=47.07; MORankSum=1.068; DD=8.81; ReadPosRankSum=0.566; SOR=0
chr13
        16003803
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=1.834; DP=6; ExcessHet=3.0103; FS=4.771; MLEAC=1; MLEAF=0.5; MQ=44.96: MORankSum=-1.834: OD=18.94: ReadPosRankSum=0.842;
chr13
                                                 113.64
                                                 37.32
                                                                  AC=2;AF=1;AN=2;DP=2;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MO=22.55;DD=18.66;SOR=0.693 GT:AD:DP:GO:PL 1/1:0.2:2:6:49.6.0
chr13
        16006318
chr13
        16007398
                                 G
                                                 37.32
                                                                  AC=2;AF=1;AN=2;DP=2;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=50;QD=18.66;SOR=0.693
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,2:2:6:49,6,0
chr13
        16008737
                                                 90.84
                                                                  AC=2;AF=1;AN=2;DP=3;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=39.03;QD=30.28;SOR=1.179 GT:AD:DP:GQ:PL 1/1:0,3:3:9:104,9,
chr13
        16009860
                                 Α
                                         G
                                                  151.64
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=-0.792; DP=7; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=48.97; MQRankSum=0.328; QD=21.66; ReadPosRankSum=1.204; SOR=
                                                 252.06
                                                                  AC=1,1;AF=0.5,0.5;AN=2;DP=7;ExcessHet=3.0103;FS=0;MLEAC=1,1;MLEAF=0.5,0.5;MQ=48.97;QD=25.36;SOR=0.941 GT:AD:DP:GQ:PL 1/2:0,4,3:7:99:269,10
chr13
        16009875
                                 CAG
chr13
        16009887
                                                 61.6
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=-0.366; DP=7; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=48.97; MQRankSum=-1.465; QD=8.8; ReadPosRankSum=0; SOR=0.446
chr13
        16011213
                                 G
                                                  166.14
                                                                  AC=2;AF=1;AN=2;DP=4;ExcessHet=3.0103;FS=0;MLEAC=2;MLEAF=1;MQ=54.78;QD=28.73;SOR=1.609
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,4:4:12:180,12,0
chr13
                                 G
        16011228
                                                 166.14
                                                                  AC=2;AF=1;AN=2;DP=4;ExcessHet=3.0103;FS=0;MLEAC=2;MLEAF=1;MQ=54.78;QD=30.97;SOR=1.609
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,4:4:12:180,12,0
chr13
        16012620
                                                 85.14
                                                                  AC=2;AF=1;AN=2;DP=4;ExcessHet=3.0103;FS=0;MLEAC=2;MLEAF=1;MQ=39.56;QD=21.29;SOR=0.693
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,4:4:12:99,12,0
                                         GA
chr13
        16015532
                                 G
                                                 67.28
                                                                  AC=2;AF=1;AN=2;DP=3;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=29.58;QD=33.64;SOR=2.303 GT:AD:DP:GQ:PL 1/1:0,2:2:6:79,6,0
                                                                  AC=2; AF=1; AN=2; DP=1; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=40; QD=27.24; SOR=1.609
chr13
        16015554
                                                  35.48
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,1:1:3:45,3,0
chr13
        16021251
                                                  70.64
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=0; DP=4; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=51.66; MQRankSum=-0.674; QD=17.66; ReadPosRankSum=0; SOR=0.69GT: A
                                                                  AC=1; AF=0.5; AN=2; BaseQRankSum=0; DP=5; ExcessHet=3.0103; FS=0; MLEAC=1; MLEAF=0.5; MQ=47.54; MQRankSum=-0.674; QD=17.66; ReadPosRankSum=0; SOR=0.69GT: A
chr13
        16021268
                                                 70.64
        16021728
                                                                  AC=2;AF=1;AN=2;DP=2;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=27;QD=18.66;SOR=0.693
chr13
                                         G
                                                 37.32
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,2:2:6:49,6,0
chr13
        16023773
                                                  119.96
                                                                  AC=2;AF=1;AN=2;DP=5;ExcessHet=3.0103;FS=0;MLEAC=2;MLEAF=1;MQ=42.69;QD=23.99;SOR=1.022
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,5:5:15:134,15,0
        16042755
                                 G
                                                  70.84
                                                                  AC=2;AF=1;AN=2;DP=3;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=31.93;QD=23.61;SOR=1.179 GT:AD:DP:GQ:PL 1/1:0,3:3:9:84,9,0
        16048306
                                                  58.32
                                                                  AC=2;AF=1;AN=2;DP=2;ExcessHet=3.0103;FS=0;MLEAC=1;MLEAF=0.5;MQ=60;QD=29.16;SOR=2.303
                                                                                                                                                             GT:AD:DP:GQ:PL 1/1:0,2:2:6:70,6,0
```

Next Step: Filter Raw (Unfiltered) VCF file


- Problem: Too many Variants
 - Non-related Human's differ by 0.1% --> 3Million SNPs (haploid genome)

- Solution: Filter Low Quality and Common Variants
 - Common Variants > 5%
 - Low Quality
 - GQ < 20, DP < 8, MQ < 40 or GATK VQSR

Filtering Criteria




VCF Visualization with IGV




Further Analysis after SNV Calling


De-novo Analysis with Trios



Rare Transmitted Analysis



Germline SNV Calling Tools



GATK4 Docker Image

Use Docker Image with GATK4 Already Installed

- Use GATK4 Docker image: https://hub.docker.com/r/broadinstitute/gatk/
- Advantage over local installation: No need to download multiple programs and worry about incompatibilities
 - Easy to use different versions of GATK4

How to use GATK4

- Use commands in same manner as Linux command line
 - Link for syntax: https://gatk.broadinstitute.org/hc/en-us/articles/360035531892
 - Includes Picard tools
- List all tools: /gatk/gatk --list

Extra Information

GATK4 Installation

Download Files

- See here for guide: https://gatk.broadinstitute.org/hc/en-us/articles/360036194592-Getting-started-with-GATK4
- Download GATK4 jar files: https://github.com/broadinstitute/gatk/releases
- Note: You will need to download a few other programs to run full SNV discovery workflows
 - BWA-MEM
 - SAM-Tools

or


Use Docker Image with GATK4 Already Installed

- Use GATK4 Docker image: https://hub.docker.com/r/broadinstitute/gatk/
- Guide for Local Docker use: https://gatk.broadinstitute.org/hc/en-us/articles/360035889991
 - Compute0 or compute1 will use slightly different syntax
- Advantage over local installation: No need to download multiple programs and worry about incompatibilities

HaplotypeCaller Modes


- HaplotypeCaller Modes
 - 1. VCF
 - For Single Sample Workflow. No need to run GenotypeGVCFs.
 - No reference confidence call for Genotypes
 - Only produces calls at variant sites
 - 2. GVCF (Default)
 - Genomic VCF File with condensed non-variant block
 - Produces calls at all sites with compressed non-variant "blocks"
 - Scales well: For joint-analysis and joint-genotyping of large cohorts
 - Best Practices Workflow
 - 3. GVCF (BP_Resolution)
 - Genomic VCF File with no compression
 - Highest Resolution, but Very large file sizes
- GVCF uses reference model to emit confidence in Genotype call
 - Updates GQ and PL annotation
 - *Intermediate File: Must be used with GATK GenotypeGVCFs in order to produce a final vcf

GVCF Format and Information



How to Install Cromwell Jar File

GO TO https://github.com/broadinstitute/cromwell/releases



Download jar file directly to HPC environment with 'wget' wget https://github.com/broadinstitute/cromwell/releases/download/53.1/cromwell-53.1.jar

Next Step Annotation

Annotate Variants with known Database Information

Annotation Tools

Databases

Ensembl Variant Effect Predictor (VEP)







