

DNA sequencing, FASTQ format, tools.

BFX Workshop Week 2

Chris Miller

Parts adapted from:

Applied Computational Genomics, Lecture 6



<https://github.com/quinlan-lab/applied-computational-genomics>

Aaron Quinlan


How to sequence a human genome: Sanger method

Key points:

- 1) sequencing by synthesis (not degradation)
- 2) primers hybridize to DNA
- 3) polymerase + dNTPS + ddNTP terminators at low concentration
- 4) 1 lane per base, visually interpret ladder



Sanger sequencing: technological advances



The next wave of DNA sequencing technologies


whiz-bang terms

- “Massively parallel” sequencing
- “High-throughput” sequencing
- “Ultra high-throughput” sequencing
- “Next generation” sequencing (NGS)
- “Second generation” sequencing


- **2005: 454 (Roche)**
- **2006: Solexa (Illumina)**
- **2007: ABI/SOLiD (Life Technologies)**
- **2010: Complete Genomics**
- **2011: Pacific Biosciences**
- **2010: Ion Torrent (Life Technologies)**
- **2015: Oxford Nanopore Technologies**

Solexa (Illumina) sequencing (2006)


- PCR amplify sample (opt.)
- Immobilize and amplify single molecules on a solid surface
- Reversible terminator sequencing with 4 color dye-labelled nucleotides



Cluster amplification by "bridge" PCR



4 different images merged



6 cycles w/ base-calling



Illumina sequencing summary

Advantages:

- Best throughput, accuracy and read length for any 2nd gen. sequencer
- Fast & robust library preparation

Disadvantages:

- Inherent limits to read length (practically, 150bp)
- Some runs are error prone



Illumina HiSeq

~3 billion paired 100bp reads
~600Gb, \$10K, 8 days
(or “rapid run” ~90Gb in 1-2 days)

Illumina X Ten


~6 billion paired 150bp reads
1.8Tb, <3 days, ~1000 / genome(\$\$)
(or “rapid run” ~90Gb in 1-2 days)

Illumina NextSeq

One human genome in **<30 hours**


Paired-end sequencing: A molecular hack to sequence longer fragments

genomic DNA



Shear to desired length (~400bp)

DNA fragments




ligate adapters, size select

sequencing library



Illumina GA2



clusters on a flow-cell




millions to billions of paired-end reads (readpairs)

~150bp ~200bp ~150bp

5' GGTGTACGAATAGTTTCCTTTTACACTCCTTGACCATCCTAGC -----//----- GGACTGAAACTTCATCTGTCTTTATAGATATGCGTGCAGCAGC 5'

-----//-----

Pacific Biosciences




Key Points:


- 1 DNA molecule and 1 polymerase in each well (zero-mode waveguide)
- 4 colors flash in real time as polymerase acts
- Methylated cytosine has distinct pattern
- No _theoretical_ limit to DNA fragment length

Caveats:

- higher error rate (1-2%)
- lower throughput : roughly 5 gigabases per run



Pacific Biosciences: long reads. Great for genome assembly



Pacific Biosciences: long reads. Great for genome assembly


Sequence and assemble a bacterial genome in 24 hours for ~ \$500
Human genome for ~\$10,000

Single haplotype assembly of the human genome from a hydatidiform mole

Karyn Meltz Steinberg,¹ Valerie A. Schneider,² Tina A. Graves-Lindsay,¹
Robert S. Fulton,¹ Richa Agarwala,² John Huddleston,^{3,4} Sergey A. Shiryev,²
Aleksandr Morgulis,² Urvashi Surti,⁵ Wesley C. Warren,¹ Deanna M. Church,⁶
Evan E. Eichler,^{3,4} and Richard K. Wilson¹

¹The Genome Institute at Washington University, St. Louis, Missouri 63108, USA; ²National Center for Biotechnology Information, National Library of Medicine, National Institutes of Health, Bethesda, Maryland 20894, USA; ³Department of Genome Sciences, University of Washington, Seattle, Washington 98195, USA; ⁴Howard Hughes Medical Institute, University of Washington, Seattle, Washington 98195, USA; ⁵Department of Pathology and Human Genetics, University of Pittsburgh, Pittsburgh, Pennsylvania 15260, USA; ⁶Personalis, Inc., Menlo Park, California 94025, USA

Oxford Nanopore Technologies




Nanopore sequencing with ONT is accurate and reliable

*load libraries onto
array of 2,048 nanopores*




2015 yield (R9 chemistry): ~500 Mbp
mean read size: 4-6 kbp




Current yield per run ("R9.4" chemistry):
~5 Gbp, 97% identity (i.e., 3% error rate)

ONT sequence length distribution





This is out of date - they're even longer now!

Nanopore sequencing is *extremely* portable



Kate Rubins sequencing DNA on the ISS

Base calling: the conversion of signal to a nucleotide sequence



Raw signal
(e.g., 454 Life Sciences)

Errors happen.
Hopefully infrequently

↓ Base calling algorithms

ACCTTCGAACGGCGGGGGGTTACAA

(Mostly) all technologies yield DNA sequences in FASTQ format

DNA



```
@seq1
ACCTTCGAACGGCGGGGGTTACAA
+
!''*(((((***+))%%%++).1***
@seq2
TGGAACCGAACGGCCCCGGTTACAT
+
!''*!!!!***+))+++++).1***
And so on...
```

The FASTQ format. Welcome to a minor hell.

A “standard” format for storing and defining sequences from next-generation sequencing technologies.

```
Sequence ID @SEQ_ID
      Sequence GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
<separator> +
Quality scores !' '* (( ( (***) ) %%%++) (%%%) .1***-+*' ')) **55CCF>>>>>CCCCCCC65
```

http://en.wikipedia.org/wiki/FASTQ_format

The FASTQ format's sequence identifier (first line of each record)

Old format

```
@HWUSI-EAS100R:6:73:941:1973#0/1
```

HWUSI-EAS100R	the unique instrument name
6	flowcell lane
73	tile number within the flowcell lane
941	'x'-coordinate of the cluster within the tile
1973	'y'-coordinate of the cluster within the tile
#0	index number for a multiplexed sample (0 for no indexing)
/1	the member of a pair, /1 or /2 (<i>paired-end or mate-pair reads only</i>)

New format

```
@EAS139:136:FC706VJ:2:2104:15343:197393 1:Y:18:ATCACG
```

EAS139	the unique instrument name
136	the run id
FC706VJ	the flowcell id
2	flowcell lane
2104	tile number within the flowcell lane
15343	'x'-coordinate of the cluster within the tile
197393	'y'-coordinate of the cluster within the tile
1	the member of a pair, 1 or 2 (<i>paired-end or mate-pair reads only</i>)
Y	Y if the read is filtered, N otherwise
18	0 when none of the control bits are on, otherwise it is an even number
ATCACG	index sequence

FASTQ quality scores: estimate of confidence in each base (sequencing technologies make errors!)

Sequence ID	@SEQ_ID
Sequence	GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
<separator>	+
Quality scores	! ' ' * (((* * * +)) % % % + +) (% % % %) . 1 * * * - + * ' ')) * * 5 5 C C F > > > > > C C C C C C C C 6 5

FASTQ quality scores: estimate of confidence in each base (sequencing technologies make errors!)

Sequence ID	@SEQ_ID
Sequence	GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
<separator>	+
Quality scores	! ' ' * (((* * * +)) % % % + +) (% % % %) . 1 * * * - + * ' ')) * * 5 5 C C F > > > > > C C C C C C C C 6 5



Qualities are based on the Phred scale and are *encoded*

$$Q = -10 \cdot \log_{10}(P_{\text{err}})$$

Phred quality score calculation

$$Q = -10 \cdot \log_{10}(P_{\text{err}})$$

Error probability (P_{err})	$\log_{10}(P_{\text{err}})$	Phred quality score
1	0	0
0.1	-1	10
0.01	-2	20
0.001	-3	30
0.0001	-4	40

FASTQ quality scores: estimate of confidence in each base (sequencing technologies make errors!)

Sequence ID	@SEQ_ID
Sequence	GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
<separator>	+
Quality scores	! ' ' * (((* * * +)) % % % + +) (% % % %) . 1 * * * - + * ' ')) * * 5 5 C C F > > > > > C C C C C C C C 6 5


Quality score encoding based on ASCII table

Dec	Hex	Char	Dec	Hex	Char	Dec	Hex	Char	Dec	Hex	Char
0	00	Null	32	20	Space	64	40	@	96	60	`
1	01	Start of heading	33	21	!	65	41	A	97	61	a
2	02	Start of text	34	22	"	66	42	B	98	62	b
3	03	End of text	35	23	#	67	43	C	99	63	c
4	04	End of transmit	36	24	\$	68	44	D	100	64	d
5	05	Enquiry	37	25	%	69	45	E	101	65	e
6	06	Acknowledge	38	26	&	70	46	F	102	66	f
7	07	Audible bell	39	27	'	71	47	G	103	67	g
8	08	Backspace	40	28	(72	48	H	104	68	h
9	09	Horizontal tab	41	29)	73	49	I	105	69	i
10	0A	Line feed	42	2A	*	74	4A	J	106	6A	j
11	0B	Vertical tab	43	2B	+	75	4B	K	107	6B	k
12	0C	Form feed	44	2C	,	76	4C	L	108	6C	l
13	0D	Carriage return	45	2D	-	77	4D	M	109	6D	m
14	0E	Shift out	46	2E	.	78	4E	N	110	6E	n
15	0F	Shift in	47	2F	/	79	4F	O	111	6F	o
16	10	Data link escape	48	30	0	80	50	P	112	70	p
17	11	Device control 1	49	31	1	81	51	Q	113	71	q
18	12	Device control 2	50	32	2	82	52	R	114	72	r
19	13	Device control 3	51	33	3	83	53	S	115	73	s
20	14	Device control 4	52	34	4	84	54	T	116	74	t
21	15	Neg. acknowledge	53	35	5	85	55	U	117	75	u
22	16	Synchronous idle	54	36	6	86	56	V	118	76	v
23	17	End trans. block	55	37	7	87	57	W	119	77	w
24	18	Cancel	56	38	8	88	58	X	120	78	x
25	19	End of medium	57	39	9	89	59	Y	121	79	y
26	1A	Substitution	58	3A	:	90	5A	Z	122	7A	z
27	1B	Escape	59	3B	;	91	5B	[123	7B	{
28	1C	File separator	60	3C	<	92	5C	\	124	7C	
29	1D	Group separator	61	3D	=	93	5D]	125	7D	}
30	1E	Record separator	62	3E	>	94	5E	^	126	7E	~
31	1F	Unit separator	63	3F	?	95	5F	_	127	7F	□

Formula for getting PHRED quality from encoded quality:

$$Q = \text{ascii}(\text{char}) - 33$$

Example:



Historically, FASTQ has had different encoding schemes for encoding PHRED quality scores. Ouch.



Current encoding:
 ! = quality 0
 J = quality 41

S - Sanger Phred+33, raw reads typically (0, 40)
X - Solexa Solexa+64, raw reads typically (-5, 40)
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 40)
with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
(Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)

FASTQE

<https://github.com/fastqe/fastqe>

```
$ fastqe example.fastq --min --max --bin
```

```
example.fastq    max (binned)
```



```
example.fastq    mean (binned)
```



```
example.fastq    min (binned)
```



Quality score binning

Table 1: Q-Score Bins for an Optimized 8-Level Mapping

Quality Score Bins	Example of Empirically Mapped Quality Scores*
N (no call)	N (no call)
2–9	6
10–19	15
20–24	22
25–29	27
30–34	33
35–39	37
≥ 40	40

FASTQ vs FASTA


>seq1

```
ACGACGAGACCTTCATCAAAAACATCATCATCCAGGACTGTATGTGGAGCGGCTTCTCGG
CCGCCGCCAAGCTCGTCTCAGAGAAGCTGGCCTCCTACCAGGCTGCGCGCAAAGACAGCG
GCAGCCCGAACCCCGCCCGCGGCCACAGCGTCTGCTCCACCTCCAGCTTGTACCTGCAGG
ATCTGAGCGCCGCGCCTCAGAGTGCATCGACCCCTCGGTGGTCTTCCCCTACCCTCTCA
ACGACAGCAGCTCGCCCAAGTCCTGCGCCTCGCAAGACTCCAGCGCCTTCTCTCCGTCCT
CGGATTCTCTGCTCTCCTCGACGGAGTCCTCCCCGCAGGGCAGCCCCGAGC
```

>seq2

```
TCCATGAGGAGACACCGCCCACCACCAGCAGCGACTCTGGTAAGCGAAGCCCGCCCAGGC
CTGTCAAAGTGGGCGGCTGGATACCTTTCCCATTTTCATTGGCAGCTTATTAAACGGGC
CACTCTTATTAGGAAGGAGAGATAGCAGATCTGGAGAGATTTGGGAGCTCATCACCTCTG
AAACCTTGGGCTTTAGCGTTTCCTCCCATCCCTTCCCCTTAGACTGCCCATGTTTGCAGC
CCCCCTCCCCGTTTGTCTCCACCCCTCAGGAATTTATTTAGGTTTTTAAACCTTCTGG
CTTATCTTACAACCTCAATCCACTTCTTCTTACCTCCCGTTAACATTTTAATTGCCCTGGG
GCGGGGTGGCAGGGAGTGTATGAATGAGGATAAGAGAGGATTGATCTCTGAGAGTGAATG
AATTGCTTCCCTCTTAACCTCCGAGAAGTGGTGGGATTTAATGAACTATCTACAAAATG
```

FASTQC: Is my sequence data any good?



LSF on compute1

- If you are not set up on compute1, do that ASAP
 - If you need help, let us know!
- Example command for an interactive job:

```
bsub -Is -M 2G -G compute-timley -n 1 -R 'select[mem>2G] rusage[mem=2G] '  
-q general-interactive -a 'docker(ubuntu:focal)' /bin/bash
```

LSF on compute1

- If you are not set up on compute1, do that ASAP
 - If you need help, let us know!
- Example command for an interactive job:

```
bsub -Is -M 2G -G compute-timley -n 1 -R 'select[mem>2G] rusage[mem=2G] '  
-q general-interactive -a 'docker(ubuntu:focal)' /bin/bash
```

Some useful UNIX commands

- `head` print the first 10 lines of a file
- `tail` print the last 10 lines of a file
getting fancy: `tail -n +2`
- `wc` count the number of characters/words/lines in a file
`wc -l` for only lines
- `less` because you don't want 3 million lines scrolling through your terminal
`q` to exit, `-S` to wrap lines (lots more useful options here)
- `grep` to search through a file (`-v` to search for lines *without* pattern)

| (pipes)

You cannot be a productive command line user until you really understand the power of pipes

```
grep TP53 genes.txt | grep "exon" | wc -l
```

This kind of construction allows you to get answers quickly!

Working with compressed data

- **tar** work with a “bundle” of data

create: **tar -cvf output.tar infile1 infile2**

extract: **tar -xvf output.tar**

- **gzip** compress a single file

create: **gzip mydata.txt** (creates mydata.txt.gz)

extract: **gunzip mydata.txt.gz** (creates mydata.txt)

Often these operations are combined

```
tar -czvf myfile.tar.gz <list of files>
```

```
tar -xzvf myfile.tar.gz
```

Sed (**s**tream **e**ditor)

Most commonly used for find and replace operations:

```
cat file.txt | sed 's/foo/bar/g' >file_fixed.txt
```

Also useful in some other contexts, like printing every 2nd line of a file:

```
sed '0~2p' file.txt
```

Schedule